

Cola®Teric 2C-LV

Amphoacetate Surfactant

INCI NAMEDisodium CCAS NUMBER68650-39-5LISTINGSUSA (TSCA);

Disodium Cocoamphodiacetate

USA (TSCA); EU (REACH); Canada (DSL); Australia (AICS); Philippines (PICCS); Korea (KECI); New Zealand (NZIoC); China (IECSC)

Cola®Teric 2C-LV is a coco imidazoline-derived dicarboxylic acid amphoteric. Cola®Teric 2C-LV is a low viscosity process modification of Cola®Teric 2C that allows for handling without the need for heat and permits cold mix formulations - and easier to empty from containers. It is more versatile and tolerant of other additives permitting greater formulation latitude.

Cola®Teric 2C-LV is one of mildest amphoteric surfactants of its class, provides a moderate lather and cleans without removing dermal lipids from the skin to give a rich, conditioning effect. It is stable over a

wide pH range and useful for institutional cleaners.

APPLICATIONS


- Shampoos
- Bath products
- · Hand and facial cleansers
- Baby wipes
- Mild personal care products
- Food preparation area cleaners
- Hospital cleaners
- Hard surface cleaners

BENEFITS

- Counter-irritants for harsh surfactants
- Excellent flash foam and lathering properties
- Effective viscosity modifiers
- Compatible with all surfactant types
- Contribute mild detergency
- Non-irritating to eyes and skin
- Couple cationics into anionic based systems
- Hard water and soap tolerant

TYPICAL PROPERTIES / STRUCTURE

Appearance	Clear Liquid	
pH, 10% aqueous	8.5	
% Solids	50.0	
NaCl, %	12.0	
Color, Gardner '98	5 Max.	
Viscosity @ 25°C	700 Max.	




Where RCO⁻ represents the fatty acids derived from coconut oil

SUGGESTED USE

Shampoos and hand soaps 5% -25%; bath products 10%-30%; body washes 5%-30%; moistened wipes 1%-4%; hospital cleaners 1%-5%; hard surface cleaners 1%-3%.

SALT PERFORMANCE TEST



FORMULATION

Ultra Mild Color Protection Shampoo

	INGREDIENT	%
1	Water	qs to 100.00
2	Suga®Nate 160 NC	36.00
3	Cola®Teric LMB	12.00
4	Cola®Lipid C	2.00
5	Mirasil [®] ADME-P	2.00
6	Pationic SSL	0.50
7	Poly Suga®Mulse D9	0.50
8	Cola®Teric 2C-LV	10.00
9	Preservative, Fragrance	qs

Procedure:

Combine Water, Suga®Nate 160NC, Cola®Lipid C, and Mirasil ADME-P in order while mixing. Begin heating to 60°C. Add Cola®Teric LMB, Pationic SSL, Poly Suga®Mulse D9, and Cola®Teric 2C-LV. Once the Pationic SSL is incorporated, cool below 50°C and add fragrance and preservative.

Typical Properties:

Appearance: Clear Liquid pH: 6.0 Viscosity: 5,000 – 7,000 cps



STORAGE / HANDLING

Cola®Teric 2C-LV should be stored in closed containers. Shipped in 55 gallon poly drums (net weight 450 lb/204.1 kg) and totes (net weight 2250 lb/1020.5 kg). Typical shelf life is 24 months from date of manufacture. Safety Data Sheets may be found at www. colonialchem.com.



225 Colonial Drive • South Pittsburg, TN 37380 www.colonialchem.com

Innovative Specialty Surfactants

Technical information contained herein is believed to be accurate. However, it is furnished without charge or obligation and is given and accepted at the recipient's sole risk. No guarantee of the accuracy of the information is made and the products discussed are sold without conditions or warranties expressed or implied. No warranties they on the guarantee that Clonaid Chemical products are manufactured to specia are expressed or implied, are the use of material is beyond the guarantee that Clonaid Actemines ustability of the product rother particular purposes. Nothing contained herein shall be considered a recommendation for any use that may infinity guora postent tight. All rights reserved to (2)(2)(2)