# Technical Data Sheet Cola<sup>®</sup>Liquid DT2



## CHEMICAL NAMETallate DIPACAS NUMBER38618-12-1EINECS NUMBER 254-043-7LISTINGSTSCA, DSL, EINECS, AICS, IECSC, ASIA-PAC, NZIoC

#### DESCRIPTION

**Cola**<sup>®</sup>Liquid DT2 is a non-ionic 2:1 diisopropanolamine (DIPA) amide of tall oil fatty acids. Colonial Chemical's Cola<sup>®</sup>Liquid products consist of completely liquid, consumer acceptable non-DEA amides that perform like traditional alkanolamide surfactants. Cola<sup>®</sup>Liquid DT2 can be used as a substitute for 2:1 DEA amides in any application where traditional alkanolamide surfactants have been used. Cola<sup>®</sup>Liquid DT2 is effective in a wide variety of household and industrial products, including:

- Synthetic, semi-synthetic metal-working compounds
- Soluable oils
- · Metalworking formulations for stamping, drawing and machining

#### **BENEFITS**

- Reserve alkalinity
- Co-emulsify mineral oils in water
- Good lubricity
- Low foam emulsions
- Stable emulsions in hard water
- Water dispersible
- Oil-in-water or water-in-oil emulsifier
- Compatibility with metalworking additives, such as corrosion inhibitors, EP agents and others
- Ease of formulating soluble oils, semi-synthetic fluids, synthetic fluids, emulsions degreasers and emulsifiable rust preventatives
- Low viscosity for ease of handling
- Contains no DEA, MEA or TEA

#### **SPECIFICATIONS**

Appearance	Clear Amber Liquid
Color, Gardner	8 Max.
pH (3% solution)	9.0 – 10.0
Acid Value	23.0 – 29.0
Alkali Value	105.0 – 131.0
% DIPA	25.0 – 31.0

#### STORAGE AND HANDLING

Cola<sup>®</sup>Liquid DT2 should be stored in closed containers. Frozen material may be thawed in a warm room, venting drums while warming. A Safety Data Sheet is available upon request.

#### PACKAGING

Cola®Liquid DT2 is shipped in 55 gal. closed-head steel drums.

### **Colonial Chemical, Inc.**

225 Colonial Drive · South Pittsburg, TN 37380 · Phone: 423-837-8800 · Fax: 423-837-3888 · www.colonialchem.com

Technical information contained herein is believed to be accurate. However, it is furnished without charge or obligation and is given and accepted at the recipient's sole risk. No guarantee of the accuracy of the information is made and the products discussed are sold without conditions or warranties expressed or implied. No warranties beyond the guarantee that Colonial Chemical products are manifadured to speas are expressed or implied, since the use of material is beyond our control. Purchasers should make their own tests and determine suitability of the product for their particular purposes. Nothing contained herein shall be considered a room room tests and the termine suitability of the product for their particular purposes. Nothing contained herein shall be considered are commendation for any use that may infinge upon patent rights. Safety information regarding this product is contained in a Material Safety Data Sheet.



1 min. 5

5 mins. 10 mins.

Shake Foam Test	Cola <sup>®</sup> Liquid DT2
1 min.	80 mm
3 min.	65 mm
5 min.	55 mm
10 min	35 mm

(1 g Test Material into 50 mL DI-water)